Manufacturer: San Huei United Co., Ltd.

Model Tested: SH 9550 Date Tested: May 28, 2020

These findings pertain to the San Huei United Co., Ltd., model SH 9550. The packaging and labeling for this product indicates that it is a NIOSH-approved product, under approval number TC-84A-3713.

Forty respirators were submitted for evaluation. However, since these are NIOSH-approved respirators, only ten samples were tested. Samples 1-5 were from the container labeled 1B37 and samples 6-10 were from the container labeled 1B39. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found here.

The maximum and minimum filter efficiency was 99.32% and 98.70%, respectively. All ten respirators measured more than 95%.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for <u>Crisis Capacity Strategies</u> (during known <u>shortages</u>).

## **Evaluation of International Respirators**



**Test:** Modified TEB-APR-STP-0059

Date Tested: May 28, 2020

Report Prepared: May 29, 2020

Manufacturer: San Huei United Co., Ltd.

Item Tested: SH 9550

**Country of Certification:** USA (42 CFR 84)


end of this report.

Pictures have been added to the

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH₂O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency	
1	85	9.5	0.68	0.68	99.32	
2	85	9.7	1.05	1.05	98.95	
3	85	9.3	0.81	0.81	99.19	
4	85	9.6	1.67	1.67	98.33	
5	85	9.8	0.74	0.74	99.26	
6	85	8.3	1.28	1.29	98.71	
7	85	9.3	0.80	0.80	99.20	
8	85	9.0	1.16	1.16	98.84	
9	85	9.3	1.26	1.27	98.73	
10	85	9.1	1.23	1.30	98.70	
	Minimum Filter Efficiency: 98.70			Maximum Filter Efficiency: 99.32		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of
  respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and
  therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.







SAN HUEI UNITEO CO., LTO. 120 Chi Pei Road, Panchiao Taipei, Taiwan Tel : (886-2) 2683-0356



#### THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:

TC -	PROTECTION 1	RESPIRATOR	CAUTIONS AND LIMITATIONS <sup>2</sup>
		SH9550	
84A-3713	N95	x	ABCJMNOP

#### 1. PROTECTION

**N95-Particulate Filter** (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.

#### 2. CAUTIONS AND LIMITATIONS

- A- Not for use in atmospheres containing less than 19.5% oxygen.
- B- Not for use in atmospheres immediately dangerous to life or health.
- C- Do not exceed maximum use concentrations established by regulatory standards.
- J- Failure to properly use and maintain this product could result in injury or death.
- M- All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N- Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O- Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P- NIOSH does not evaluate respirators for use as surgical masks.



#### A WARNING

NIOSH APPROVED:N95 95% efficient against particulate aerosols free of oil.

This respirator helps protect against certain particulate contaminants but does not alliminate exposure to or the risk of control that disease or infection. Misuse many materials in sickness or death.



### **A** ADVERTENCIA

APPROBADO POR NIOSH:N95 Ofrece 95% de eficiencia contra aerosol de partículas libre de aceite.

Este respirador ayuda e proteger contra ciertos polvos y neblinas pero no elimina la exposición o el riesgo de contraer infecciones o enfermedas. EL uso incorrecto de este producto puede ser causa de enfermedad o muerte.















