NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

Manufacturer: Inovel, LLC Model Tested: 3002N95-M Date Tested: May 27, 2020 Report Prepared: May 28, 2020

These findings pertain to the Inovel, LLC, model 3002N95-M (TC-84A-4102) submitted for testing, and may not be applicable to other stockpile facilities and/or under different environmental storage conditions. The maximum and minimum filter efficiency was 99.64% and 98.88%, respectively. All twenty respirators measured more than 95% efficiency. No shelf life was designated for this model by the previous NIOSH approval holder.

NIOSH regulation sets the minimum quality and performance requirements for the approval of respirators (42 CFR 84). NIOSH does not have requirements for shelf life or storage conditions for particulate-only APRs. The approval holder (i.e. the entity that is granted the approval from NIOSH) is responsible for understanding how their products' design or performance may be affected by various use or storage conditions and must provide instruction for establishing the proper use, storage, and maintenance procedures for their approved products, which may include designating a shelf life. FFR or particulate filter packaging (such as the box) often includes NIOSH-approved user instructions, label information, and recommendations on shelf life. Additionally, some approval holders also disseminate recommendations related to storage and shelf life through resources such as user and web notices. The respirators tested in this study were generally not designed for long-term storage.

Based on research conducted by NIOSH and this limited testing, NIOSH does not have enough information to definitively know the level of protection that may be provided by respirators that 1) are stored for prolonged periods of times; 2) are stored under various storage conditions; or 3) have exceeded the approval holder's designated shelf life. Users of respirators that have exceeded the designated shelf life should be forewarned to avoid a false sense of confidence; these devices may not provide the same level of protection as those that have not exceeded the designated shelf life. We recommend contacting the approval holder(s) of the respirators in the stockpile with specific questions regarding the use of product beyond the manufacturer- designated shelf life.

The results provided in this letter are specific to the subset of NIOSH-approved N95s, past their designated shelf life, that were provided to NPPTL for evaluation.

These results will be added to the CDC guidance for <u>Stockpiled N95 Filtering Facepiece Respirators</u> <u>Beyond the Manufacturer-Designated Shelf Life</u>.

Evaluation of Stockpiled and Beyond Manufacturer-Designated Shelf Life N95s


Test: TEB-APR-STP-0059

Manufacturer: Inovel, LLC Item Tested: 3002N95-M

Expiration Date: None Provided **Manufacture Date:** None Provided

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	11.5	1.09	1.09	98.91
2	85	17.2	0.65	0.65	99.35
3	85	14.3	0.63	0.63	99.37
4	85	11.3	0.93	0.93	99.07
5	85	11.9	0.84	0.84	99.16
6	85	16.7	0.72	0.72	99.28
7	85	15.5	0.75	0.75	99.25
8	85	12.2	0.79	0.79	99.21
9	85	17.8	0.36	0.36	99.64
10	85	16.7	0.55	0.55	99.45
11	85	11.7	0.90	0.90	99.10
12	85	15.9	1.12	1.12	98.88
13	85	14.0	0.71	0.71	99.29
14	85	11.4	0.95	0.95	99.05
15	85	15.8	1.08	1.10	98.90
16	85	11.9	0.77	0.77	99.23
17	85	16.7	0.51	0.51	99.49
18	85	13.9	0.65	0.65	99.35
19	85	17.0	0.57	0.57	99.43
20	85	11.8	1.04	1.04	98.96
Minimum Filter Efficiency: 98.88			Maximum Filter Efficiency: 99.64		

DESCRIPTION:

The Inovel 3000 Series Health Care N95 Particulate Respirator and Surgical Mask is designed to help provide respiratory protection for the wearer. This product has been tested and certified by NIOSH as an N95 respirator and as having a filter efficiency level of 95% or greater against particulate aerosols free of oil. It is fluid resistant, disposable and may be worn in surgery or throughout the hospital.

This product contains no components made from natural rubber latex.

INTENDED USE:

The various models of Inovel Type N95 Healthcare Particulate Respirators and Surgical Masks are intended to be worn by operating room personnel during surgical procedures to protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.

The devices also meet CDC Guidelines for TB Exposure Control.

CONTRAINDICATIONS

- 1. Not for use with beards, or other facial hair that prevents direct contact between the face and sealing surface of the respirator.
- 2. Eyewear must not prevent direct contact between the face and sealing surface of the respirator.
- 3. Not to be used on children.
- 1 Tested in accordance with NIOSH 42 CFR 84 2 Passed ASTM F 1862 @ 160mm Hg

LIMITED WARRANTY – IMPORTANT NOTICE TO PURCHASER

This limited warranty is made in lieu of the warranties of merchantability. fitness for particular purposes and all other warranties, express or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Inovel will be met by products sold. Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Inovel products proved to be defectively manufactured. Except as provided above, Inovel shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

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
U.S. Patents #D285,374, #4,850,347 and #4,384,577 Trademarks and Patents of Moldex-Metric, Inc., Culver City, CA 90232 USA 3000-802 REV B 4/06 Made in U.S.A.

FITTING

- If you cannot obtain a proper fit do not enter the contaminated area.
 OSHA regulation 29 CFR 1910.134(f) requires that the user must be fit tested before being assigned a respirator. Use Inovel Qualitative Fit Test Kit or an OSHA approved Quantitative fit test such as TSI's PortaCount.
- Users must follow instructions each time respirator is worn. It's the wearer's responsibility to obtain a proper fit.
- To choose the correct 3000 N95
 Series size/ shape you must fit test and user seal check.
- When changing from one Inovel model or size to another or when changing from another manufacturer's product to Inovel, you must fit test with the new respirator.
- Not for use with beards, or other facial hair that prevents direct contact between the face and sealing surface of the respirator.
- Eyewear must not prevent direct contact between the face and sealing surface of the respirator.

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3000 SERIES

- 1. Hold respirator in hand with molded nose contour (narrow end) at finger tips, allowing headstraps to fall below hand.
- 2. Place respirator under chin with molded nose contour (narrow end) up. Cushion, if included, must not be creased inside respirator. Raise top strap to top back of head. Pull shorter bottom strap over head, below ears, to around neck. Do not wear with only one strap because it may affect fit.
- 3. Adjust respirator for comfortable fit.
 4. Each time user enters contaminated area seal should be fit checked. Cover front of respirator by cupping both hands. INHALE SHARPLY. A negative pressure should be felt inside respirator. If any leakage is detected at respirator edges, adjust straps by pulling back along the sides and/or reposition respirator. Repeat until sealed properly, otherwise see your supervisor. Entry into a contaminated area with an improper fit may result in sickness or death

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A WARNING TO USER:

Non-latex straps are used in the headbands of these respirators. However, individuals highly sensitive to natural rubber latex may have an allergic reaction.

This product does not eliminate the risk of contracting any disease or infection. Change immediately if the respirator comes in contact with blood or fluids, including body fluids. OSHA has not set a permissible exposure level for airborne biohazards.

This product has not been sold with warning or use instructions for personnel involved in industrial or related situations.

If you are considering such uses, first contact the Inovel Technical Services pepartment at 866-546-6835, ext. 523.

RESTRICTIONS:

- Before use, a written respiratory protection program in accordance with 29 CFR 1910.134 must be implemented.
- 2. This respirator does not supply oxygen and must not be used in atmospheres containing less than 19.5% oxygen.
- Prior to each use, inspect filter media and strap attachment area for tears and damage. Staple or strap attachment perforations do not affect NIOSH approval.
- If respirator is damaged, distorted, a proper fit cannot be obtained or breathing becomes difficult, leave contaminated area as soon as possible and replace respirator.
- 5. If the respirator comes in contact with blood or fluids, including body fluids, leave contaminated area as soon as possible and discard and replace the respirator,
- 6. Do not alter, modify or abuse this respirator.
- 7 Store unused respirators in box in a clean, dry, non-contaminated area.
- Dispose of respirator according to your employer's policy and local regulations.
- Do not reuse or store for reuse or harig around neck unless your employer specifically authorizes reuse.
- 10. When used for surgical procedures, discard after every use.
- 11. Use respirator before the "use by" expiration date printed on box.
- For technical assistance call Inovel Technical Service Department, 866-546-6835, ext. 523.

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